

TAMPER EVIDENT DOPING CONTROL SOLUTIONS FOR CLEANER SPORTS

INTEGRITY • TAMPER EVIDENCE • HERITAGE

VERSAPAK HAVE BEEN MANUFACTURING HIGHLY SECURE ANTI-DOPING EQUIPMENT IN THE UK SINCE 1987, OUR EQUIPMENT IS USED EXCLUSIVELY IN THE WORLDWIDE FIGHT AGAINST THE USE OF DRUGS IN SPORT.

Recent years have shown us that the fight for clean sport has become ever more complex. Doping is widespread and threatens fair competition. The sophistication used to circumvent current Doping Control standards is increasing, to the point the best cheats can win, not the best athlete! Therefore making ongoing improvements in the sample collection process critical, and Versapak's mission.

At Versapak Doping Control, our aim is to ensure that our Anti-Doping testing equipment meets the needs of those working towards fair competition in sport. We partner with sporting events, federations, national and international Anti-Doping organisations and WADA-accredited laboratories worldwide.

Our message throughout this product catalogue is of our core values "Integrity – Tamper Evidence – Heritage." We strive to deter those looking to undermine the spirit of clean sports as there is no place for cheating, in any form, true recognition is one that is honestly earned. Versapak manufacture tried and tested tamper evident products in the UK that are used in countless other industries and our expertise in tamper evidence reaches back to 1973 when Versapak was first established.

We respect the athletes' right to complete assurance in the Anti-Doping process, and understand that athletes require total confidence in the process. That's why, unlike our competitors, we send our sample kits to internationally accredited packaging organisations, who specialise in testing secure and tamper evident packaging products, to evaluate their integrity and ensure they are fit-for-purpose. Their findings have concluded that Versapak's Sample Collection Kit "cannot be accessed without showing evidence of tampering."

I am confident you will find our products will be of interest to your organisation and that our products and values align with your views of protecting clean sport. My team and I look forward to supporting your Doping Control programme with Versapak.



CAROLINE ATKINSON Group Managing Director



INTEGRITY

and not a fraudulent copy.

and visible signs of tampering.

HERITAGE

Versapak anti-doping equipment is specifically designed to comply with the World Anti-Doping Agency (WADA) requirements on sample collection equipment, as defined in the WADA International Standards Testing & Investigations (ISTI) of March 2019.

SINCE OUR BEGINNING, WE HAVE VALUED CLEAN SPORT, FAIR PLAY, TAMPER-EVIDENCE AND INTEGRITY - AND THE SAME STILL BEATS STRONG IN OUR HEART TODAY.



Our blood and urine kits are manufactured with a uniquely identifiable DNA marker. The same technology is used in bank notes within the banking industry to determine authenticity. It is detectable with a scanner to offer reassurance to users and the industry that our kits are a genuinely manufactured Versapak Doping Control Kit

TAMPER EVIDENCE

Our kits are designed to reveal any interference with the contents. Unlike competitors, Versapak considers this fact above-all-else to be the defining aspect of protecting the integrity of the athletes sample. We appointed an independent world renowned secure packaging organisation who thoroughly tested the tamper evident credentials of our product. They concluded that the kit could not be accessed without showing clear



Our expertise in Tamper Evident Security products dates back to 1973, so it's fair to say we have a wealth of experience in the field of security. We worked alongside organisers of the 1987 Rugby World Cup with their first ever Doping Control tests. We are proud of our heritage and experience in Doping Control.





PERSONALISATION

Whether you're looking for added security or some event advertising, why not personalise your blood or urine kits with your own information or logo? By sending us your high resolution image, we can laser print onto the canisters so that it can never be scratched or removed.



WIDE OPENING AIDS USABILITY

A wide neck opening makes filling the sample bottle to the correct level much less problematic than others, particularly when doing so with one of our easy-pour collection cups.

MEDICAL GRADE 'PETG' BOTTLE

Better in many ways, the medical grade material is strong and robust and is particularly suitable for biological samples. Tests show the bottle can be frozen over 3 freeze/thaw cycles, without succumbing to thermal shock, thanks to its superior structural integrity over glass.



UNBEATABLE SECURITY

Our patented once-only locking mechanism ensures that a Versapak Urine and Blood canister cannot be reopened without showing clear evidence of tampering.

DOUBLE SKINNED PROTECTION

The two layers of the bottle and canister makes it impossible to gain access to the sample without

ALL OUR EXPERTISE HAS GONE INTO THE PERFECTION OF THIS TAMPER EVIDENT **PRODUCT THAT ADHERES TO** WADA, INADO AND SPORTS **INDUSTRY STANDARDS**





DNA SCANNER

Our blood and urine kits both include a unique DNA marker that works with our security scanner, offering reassurance to laboratories that the kit is a genuine Versapak product.



URINE COLLECTION VESSEL

A simple and effective product for collecting the urine sample. Collect and seal in a leak proof container until ready to pour into the sample bottles. Features fill level markers, up to 125ml, making it easier to decant the correct amount into the A and B bottles. Also available as a partial urine storage kit, contact Versapak for more details.

Item no: UCV

ALL VERSAPAK URINE SAMPLE KITS ARE MANUFACTURED TO MEET WADA ISTI GUIDELINES. OUR NEWEST INNOVATION IS THE STRONGER, MEDICAL GRADE, PETG BOTTLE THAT CAN FREEZE TO -80°C OVER THREE FREEZE THAW CYCLES

URINE SAMPLE COLLECTION



The Versapak secure holdall with tamper evident security sealing is used to carry and protect urine samples and maintain the internal temperature for up to 6 hours which can be monitored with a temperature monitor. The medium holdall can carry up to 10 A/B Urine Kits in their boxes. Also available in small and large sizes.

Item no: TRANSPORTBAG1/2/3 NUMBEREDT2

A/B URINE

SAMPLE KIT

Our newest innovation is the super

strong, medical grade, PETG bottle that

can freeze to -80°C and features a DNA

marker, double skin, and the patented

winner

'once-only' locking mechanism with

unique serial identification numbers

for tamper evident security.

reddot award 2019

Item no: DTC3



PARTIAL URINE STORAGE KIT

Our partial sample collection kit is for when only a partial sample is given by the athlete and a full sample is still needed. Its wide neck makes it easy to pour, it has a tamper evident label and bag for secure storage between sample receipts.

Item no: BPK3



Unlike other kits and equipment, all Versapak products comply with WADA ISTI 2019

BLOOD COLLECTION ACCESSORIES

Our Collection Accessories contain the needles and blood sample tubes necessary for blood sample collection.

ltem no.	BD Vacutainer®	BD Security Needle 21G	BD Butterfly Needle 21G
Collection Kit 1A	2 x 5ml		
Collection Kit 1B	2 x 5ml		
Collection Kit 2A	2 x 3ml		
Collection Kit 2B	2 x 3ml		
Collection Kit 3A	2 x 3ml		
Collection Kit 3B	2 x 3ml		
Accessory packs for blood	Description		
Accessory pack ACC1	Provides sterile disinfectant pads, cotton swabs, gloves & plaster		



rr CO

INSULATED SECURE HOLDALLS FOR TRANSPORTING KITS

The Versapak secure holdall is tamper evident with security sealing. It is used to carry and protect blood samples and maintain the internal temperature for up to 6 hours which can be monitored with a temperature monitor. The small holdall can carry up to 14 A/B Blood kits in their boxes. Also available in medium and large sizes.

Item no: TRANSPORTBAG1/2/3 NUMBEREDT2

A/B BLOOD SAMPLE KIT

Featuring a DNA marker to guarantee a genuine Versapak product, double skin security, our patented 'once-only' locking mechanism and unique serial identification numbers for tamper evident security.

Item no: BLDAB

versapak

Blood Sample Transport Kit






BIOLOGICAL PASSPORT KIT

The Blood transport kit for the biological passport is used for the athlete tests according to the section 9 and Annex K4 of the International Standards of Testings' and Investigations and uses the Versapak patented tamper evident security locking system.

Item no: BLDBIO



BLOOD SAMPLE COLLECTION



TEMPERATURE MONITORS

Versapak temperature loggers from Tempmate[®] are used for monitoring temperature sensitive products throughout the entire cold chain. Available in both single-use and multi-use data formats.



EASE-OF-USE TO THE CORE

The Urine Opening tool is designed for ease-of-use, while ensuring integrity and validity when processing the samples. Once placed in its holster the DNA scanner can quickly determine genuine Versapak kits are being opened. The two handles make it effortless to turn and destroy the lid of the canister and access the sample bottle.



OPENING TOOL FOR URINE

The Versapak Urine Kit opening tool is made easy with two, lever handles for twisting and rotating the canister lid from its base. The tool is designed to carefully open the canister without affecting the bottle inside, while ensuring complete breakage so that the canister cannot be resealed covertly. It also features a holster for the DNA scanner for a streamlined opening process.

RESEAL CANISTER

Useful when needing to reseal both samples for further testing within a laboratory. Also featuring a DNA marker, patented 'once-only' locking and unique serial identification for tamper evident security.





SUPPORTING THE ANTI-DOPING **INDUSTRY LABORATORIES** WITH TOOLS TO OPEN SAMPLE KITS AND A PROCESS TO CHECK

FOR EVIDENCE OF TAMPERING

LABORATORY

OPENING TOOL FOR BLOOD

The Versapak Blood Kit opening tool is an effective design that allows for a simple corkscrew opening of the blood kit. The operation will crack the lid, protecting the vial inside so the canister cannot be covertly resealed.

DNA SCANNER

The DNA scanner is used by laboratories to confirm the authenticity of the doping control kits. With one scan, the invisible DNA marker is detected to confirm a genuine Versapak canister.

Versapak Doping Control develop sample testing kits for various animal sporting disciplines.

Our expertise in tamper evidence for over 45 years has allowed us to spearhead a new initiative to protect the integrity of animal sports events. We provide sample testing kits for Blood and Urine, as well as Collection Accessories and Transport Bags for Horse, Camel and Dog Racing events.

FOR MORE INFORMATION OR TO ASK FOR A QUOTE, PLEASE CONTACT US.

> ANIMAL SAMPLE TESTING EQUIPMENT FOR VARIOUS ANIMAL SPORTING DISCIPLINES TO KEEP ALL SPORTS CLEAN

DOPING CONTROLS FOR ANIMALS

Whether you're an Anti-Doping organisation, or a laboratory tasked with conducting the sample tests, we provide a knowledge base to support the fight against doping-in-sport.

Our content includes crucial information for laboratories to safely open Versapak Sample Kits and how to detect evidence of tampering.

Visit: **versapak-anti-doping.com** Register to gain access to restricted content.



WEBSITE – REGISTER WITH VERSAPAK DOPING CONTROL TO ACCESS ANTI-DOPING ORGANISATION AND LABORATORY CONTENT



Versapak Doping Control support WADA's Anti-Doping strategy and the Anti-Doping code of conduct.

Versapak follows and has contributed to the rules set-out in the International Standards of Testing & Investigations and Olympia Decree as set out at the International Olympia Symposium.

We have also worked with many national anti-doping organisations and laboratories such as UK Anti-Doping, Kings College, Cologne Laboratory and we have engaged with the athlete movement to pave the way to a brighter future for all stakeholders of Sports Integrity.

Follow our journey on Twitter, LinkedIn, YouTube & Instagram.



VERSAPAK IN THE MEDIA – WE CONNECT WITH ATHLETES AND ANTI-DOPING ORGANISATIONS TO DISCUSS THE ISSUES THAT MATTER AS WELL AS ACTING AS A CRITICAL VOICE OF AUTHORITY ON DOPING CONTROL TOPICS

SUPPORT

Versapak Doping Control Solutions

10

1

4 Veridion Way • Erith • Kent • DA18 4AL **T:** +44 (0)20 8333 5320 E: info@versapak-anti-doping.com

versapak-anti-doping.com