



# Athlete Information

*Not seen the Versapak Urine Collection Kit before? See how it works!*

1




Fill collection vessel to required amount.

2




Open red bottle "A" and fill to required amount.  
Repeat for blue bottle "B".

3




Ensure bottle lids are tightly sealed.  
Place bottles into corresponding security canisters.

4



Push down lid firmly until you hear clicks.  
Locking legs will now be locked in place.



More information can be found on our website, register to download instructions or why not watch our videos on our Twitter & YouTube channel



Scan for 'how to use' video





## Focus on your Finish

Versapak Doping Control kits are derived from our core values:

Integrity

Tamper Evidence

Heritage

We are convinced that our products are the best solution for combatting doping in sport. A promise that we will uphold standards for our champions so that you don't need to.




Trusted UK manufacturer since 1987



Supplying globally to over 30 countries



- ✓ Award winning kit
- ✓ Unbeatable security
- ✓ DNA Marker Technology
- ✓ Medical Grade 'PETG' material
- ✓ Double Skin Protection
- ✓ Wide Opening Neck



Versapak work with athletes to champion clean sport. Watch our #athleteintegrity interviews.



**Scan to watch our athlete interviews**



Contact  
info@versapak-anti-doping.com  
T: +44 (0)20 8333 5320

Versapak Doping Control  
Unit 4, Veridion Way  
Erith, Kent, DA18 4AL  
www.versapak-anti-doping.com